Bodyclock

Classic with Radio 150

OPERATING INSTRUCTIONS

Bodyclock Classic with Radio 150

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Bodyclock Classic mit Radio 150

BRUKSANVISNING

Bodyclock Classic med Radio 150

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OPERATING INSTRUCTIONS

Bodyclock Classic with Radio 150

- 1 Safety
 - 2 Getting started
 - 3 The controls
 - 4 Setting the time
 - 5 Setting the sunrise (to wake you up)
 - 6 Back-up audible wake-up signal
(radio or beeper)
 - 7 Using the radio
 - 8 Setting the sunset (go-to-sleep program)
 - 9 Using the Bodyclock as a reading light
 - 10 Problem solving
-

INTRODUCTION

Congratulations on buying a Bodyclock Classic with Radio. Bodyclocks work by using your body's natural response to sunrise and sunset. With proper use, it will help regulate your circadian rhythm (sleep/wake pattern), leaving you feeling more energetic. Thousands of people have benefited from its positive effects, and it is especially useful in alleviating the symptoms of the Winter Blues or SAD (Seasonal Affective Disorder). In clinical trials, dawn-simulators like the Bodyclock Classic with Radio have been shown to improve mood, energy, productivity and the quality of sleep and awakening.*

Please read these instructions carefully to get the most out of your Bodyclock Classic.

*Acta Psychiatr Scand 1993 Jul; 88 (1):67-71

1 Safety

- Keep the globe over the bulb, unless the bulb is cold and the Bodyclock unplugged.
- Don't cover the globe as the slots in it help lose heat from the bulb. The globe will get warm so keep it clear of anything that might deteriorate in contact with heat.
- Always use a E14 (SES) candle bulb (60W max).
- Keep the Bodyclock away from water and damp. For indoor use only.
- Don't open the casing for any reason. There are no user-serviceable parts inside.

2 Getting started

- Remove the 'globe', the cover over the bulb, with a gentle upward tug.
- Screw in the bulb firmly and replace the globe.
- Find a good position for the Bodyclock close to the bed-head.
- Plug the Bodyclock into the mains.
- When replacing the bulb, always unplug the Bodyclock. Spare sunshine neodymium bulbs are available from the Lumie Careline. You may need to set the time again after you plug it back in.


3 The controls

- **L** turns the dial illumination on and off.
- **-** turns the light down or off.
- **+** turns the light up.
- **A** starts the sunset and sets the wake-up sequence. It also stops the wake-up sequence in the morning.

The slide switch at the side has four positions:

- **RADIO ON** turns the radio on
 - **WAKE UP to RADIO** makes the radio come on after the sunrise, as a back-up reminder that it's time to get up.
 - **WAKE UP to BEEP** makes the beeper come on after the sunrise, as a back-up reminder that it's time to get up.
 - **WAKE UP to LIGHT** turns the radio and beeper off so there is no sound after the sunrise.
-
- **AM/FM** switch is used to select the waveband.
 - Tuning dial tunes the radio.
 - Volume dial controls the volume of the radio.

4 Setting time

- Pull out the small knob at the back of the Bodyclock to set the time. This knob is also used to set the wake-up time.
- **IMPORTANT:** make sure you turn the knob in the direction of the arrow next to it, or you may damage the clock mechanism.
- When the clock shows the correct time, push the knob back in.

5 Setting the sunrise (to wake you up)

- Use the small knob at the back of the Bodyclock to set the orange hand to the desired wake-up time.
- **IMPORTANT:** make sure you turn the knob in the direction of the arrow next to it, or you may damage the clock mechanism.
- **NOTE:** the wake-up time is the time at which the light reaches full brightness. If, for example, you wanted to wake at 7:00, point the orange wake-up hand to 7:00; at 6:30 the sunrise will start and the light will come on with a very faint glow, getting slowly brighter until it reaches full brightness at 7:00.
- When **A** is up the wake-up function is ON.
- When **A** is down the wake-up function is OFF.

6 Back-up audible wake-up signal (radio or beeper)

- Use the slide switch at the side to set an audible signal to come on after the sunrise has reached full brightness.
- To have the radio come on at wake-up time, after the sunrise, set the switch to **RADIO**.
- To have the beeper come on at wake-up time, after the sunrise, set the switch to **BEEP**.

- If you don't want an audible back-up signal, set the switch to **LIGHT**.
- To turn off the radio or beeper after it has sounded, press **A**.
- If the radio or beeper is not turned off manually, after a while they will automatically shut off, together with the light.

7 Using the radio

- You can listen to the radio during the day or at night before you are ready to go to sleep.
- To turn the radio on, set the slide switch to **RADIO**.
- When you want to turn the radio off, move the switch to another setting.
- Use the waveband switch to select **AM** or **FM**.
- The AM aerial is built in, so you may have to turn the Bodyclock around slightly to find the best reception.
- The FM aerial is the wire coming from the back of the Bodyclock. Try moving the wire around to find the best reception.
- The numbers printed on the inside of the tuning dial show the AM frequency; the outer scale shows the FM frequency.

8 Setting the sunset (go-to-sleep program)

- Press **+** until the light is as bright as you want it.
- Press **A** so the button is up.
- The Bodyclock will fade down over 30 minutes to darkness. In the morning, the wake-up sequence will activate at the set time.

9 Using the Bodyclock as a reading light

- Press + until the light is as bright as you want it.
- Press – to dim or switch off the light completely.

10 Problem solving

- **The light will not come on but the clock still works.**

Check to make sure the bulb is fully screwed into the socket. If this fails, try replacing the bulb.

- **Neither the clock nor the light work, or the Bodyclock is not doing what I expect it to.**

Unplug the Bodyclock, wait several seconds (count to 10) then plug it back in. Sometimes the microprocessor picks up static electricity and unplugging it removes the static. If this fails, check the fuse in UK plugs.

If you have any other problems with your Bodyclock, please call the Lumie Careline (0845 658 9292) before sending the unit back for repair.

Sleep advice

A good night's sleep is key to your mental and physical health. Achieve it by:

- Using a Bodyclock
- Having curtains that block out light
- Avoiding eating late at night
- Drinking a warm milky drink

Bodyclock Classic Guarantee

This unit is guaranteed against any manufacturing defect for 24 months from the date of purchase. This is in addition to your statutory rights. If you are returning your Bodyclock Classic for service, it must be in its original packaging including internal pieces. Please call us for further advice.

Contact us

Outside In is Europe's leading specialist in light therapy. Our products are derived from the latest scientific studies and we are specialists in sleep science, circadian rhythms and SAD (Seasonal Affective Disorder).

If you have any suggestions that you feel could help us improve our products or service, or if there are problems you think we should know about, please tell us.

- **Lumie Careline 0845 658 9292**
- **www.lumie.com**

English

Technical specification

Specification subject to change without notice.

Lumie Bodyclock Classic with Radio model NDE/NDU

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© Outside In November 2005
Lumie is a trademark of:
Outside In (Cambridge) Ltd